

Sample requirements

4.5 ml peripheral blood in EDTA

Incomplete forms will result in delay of sample analysis, and possible rejection of samples.

Test required

BCR-ABL Monitoring (must be received in laboratory within 48 hours of draw)

CML Mutation Analysis (must be received in laboratory within 48 hours of draw)

TPMT Genotyping

OTHER

Patient information

Surname

Forename

DOB

Male / Female

NHS No/Ref No

Postcode

Referring clinician

Name (print clearly)

Hospital

Contact Tel

Email

Sample information

Date sample taken

Sample ref

Copy report to be issued to

Name

Hospital

Email

Clinical Information

Clinical Diagnosis

Other Relevant details

Diagnosis of CML

Y/N (if appropriate)

Year of Diagnosis

Patient Taking

Imatinib / second generation TK Inhibitor

Samples should be sent to:

Sample Reception, Northern Genetics Service, Biomedicine East Wing, International Centre for Life, Newcastle upon Tyne, NE1 3BZ

Test result enquiries: 0191 241 8804 / 241 8787 (Opening hours : 8:30am – 5.00pm, Monday to Friday)

Shipping Instructions: All samples should be appropriately packaged and labelled and sent via Courier or Royal Mail 1st Class Special Delivery

For Laboratory use only

Sample ref:

Sample type