## Haemato-Oncology Testing



Sample requirements								
4.5 ml peripheral blood in EDTA								
Incomplete forms will result in delay of sample analysis, and possible rejection of samples.								
Test required								
	BCR-ABL Monitoring (must be received to NewGene within 48 hours of draw)							
	CML Mutation Analysis (must be received to NewGene within 48 hours of draw)							
	TPMT Genotyping							
	OTHER							
Patient inf	atient information		Referring clinician					
Surname			Name (print clearly)					
Forename			Hospital					
DOB		Male / Female						
NHS No/Ref No			Contact Tel					
Postcode			Email					
Sample information			Copy report to be issued to					
Date sample taken			Name					
Sample ref			Hospital					

Clinical Information				
Clinical Diagnosis				
Other Relevant details				
Diagnosis of CML	Y/N (if appropriate)			
Year of Diagnosis				
Patient Taking	Imatinib / second generation TK Inhibitor			

Email

Samples should be sent to:

NewGene, Sample Reception, Institute of Genetic Medicine, Biomedicine East Wing, International Centre for Life, Newcastle upon Tyne, NE1 3BZ.

Test result enquiries: 0191 242 1923 / Email: tnu-tr.newgene@nhs.net (Opening hours: 8:30am – 5.00pm, Monday to Friday)

Shipping Instructions: All samples should be appropriately packaged and labelled and sent via Courier or Royal Mail 1st Class Special Delivery



For NewGene use only								
Sample ref:		Sample type						